

# rainbow Acoustic Monitoring® (RAM™) Sensor and Cable Compatibility Requirements

## Respiratory Acoustic Sensors (RAS) and RAM Dual Cables

- > Masimo rainbow SET™ technology V7.6 or higher
- > RRa® parameter loaded on the device
- > Masimo RAM Dual Cable
- > Masimo rainbow® RAS sensor
- > SpO2 monitoring with a Masimo SET® or Masimo rainbow SET sensor

For professional use. See instructions for use for full prescribing information including indications, contraindications, warnings and precautions. Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

© 2016 Masimo Corporation. All rights reserved.

**Masimo U.S.**  
Tel: 1 877 4 Masimo  
info-america@masimo.com

**Masimo International**  
Tel: +41 32 720 1111  
info-international@masimo.com

